

The Detection of Pregnancy Associated Glycoproteins (PAGs) in Routine Milk Recording Samples as an Indicator of Pregnancy in Dairy Cattle


Todd Byrem, Antel BioSystems

Kathy Velek & Hannah Pearse, IDEXX Inc



RUMINANTS
WIEDERKÄUER
RUMINANTS
RUMIANTES
反刍动物
反芻動物

The Importance of Reproductive Efficiency in Dairy Herds



<http://www.extension.org/pages/11007/pregnant-vs-open:-getting-cows-pregnant-and-the-money-it-makes>

Test With Confidence™




Pregnancy Diagnosis

- Rectal Palapation/Ultrasound
 - Advantages
 - Accurate, real-time, pregnancy plus....
 - Disadvantages
 - Labor and technically intensive, disruptive to pregnancy
- Blood Chemical Assays
 - Progesterone / ECF/ estrone sulphate
 - Pregnancy Associated Glycoproteins
 - Advantages
 - Accurate, simple
 - Disadvantages
 - Labor intensive
- Milk Chemical Assays
 - Progesterone / estrone sulphate

IDEXX Milk Pregnancy Test

A New Opportunity for Laboratory-Based Pregnancy Testing

- **Accurate determination of pregnancy status** - High levels of sensitivity and specificity, from 60 days after calving. Day of detection for commercial kit to be confirmed
- **Trusted, timely results** - Obtain results in less than 3.5 hours using proven IDEXX ELISA technology
- **Expanded testing options** - Test for pregnancy from routine DHI milk samples
- **Improved reproductive performance** - Through earlier identification and rebreeding of open cows, shorter calving intervals and increased milk production



Test With Confidence™

IDEXX

IDEXX Milk Pregnancy Test

Kit Overview




- ELISA strip well format
- DHI whole milk samples
- ~3.5 hours to results
- Ready to use reagents
- Standard laboratory equipment
- Detects pregnancy associated glycoproteins (PAGs)

IDEXX Milk Pregnancy Test

Kit Overview

Pregnancy Associated Glycoproteins (PAGs)

- Target antigen for the IDEXX Milk Pregnancy ELISA
- Placenta-specific expression
 - Expressed in maternal & embryonic regions of placenta
- Subgroup of Aspartic Protease family
 - 22+ bovine transcribed genes identified
- Temporally expressed
 - Variable gene expression at different stages of pregnancy




Guruprasad, K et al. 1996 Protein Engineering. 9:849

Test With Confidence™

IDEXX

PAG Levels in Serum

PAGs are present throughout pregnancy




PAGs present throughout pregnancy. Early detection (from 28 days post AI).
PAGs decline post calving: test returns to negative from 60 days post calving

Test With Confidence™



IDEXX Milk Pregnancy Test

ELISA Protocol




Assay Format:

- Sandwich ELISA
- Biotin & SA-HRPO detection
- Color development = pregnant

Temporal Study

PAG Levels Throughout Pregnancy

- High degree of variability in PAG levels detected from different cows
- PAGs are detectable early in pregnancy and throughout gestation
- Assay response declines slightly between day 40 through day 90
- Very strong signal in late gestation through calving




5-13 cows sampled at each time point; average response and standard deviation shown on graph

IDEXX Milk Pregnancy Test

PAG Levels Post-Calving

134 cows sampled after calving, and prior to next breeding;
Assay response for individual cows is shown on the graph



- PAG levels decline rapidly after calving
- Specificity of **100%** by **60 days** (8.5 weeks) after calving
- No interference when testing for the next pregnancy

Test With Confidence™



IDEXX Milk Pregnancy Test

Results – Herd A

New England Dairy : 192 cows sampled

All cows > 60 DIM and not bred or >40 days post insemination

		Ultrasound or Palpation	
		Pregnant	Open
IDEXX Milk ELISA	Pregnant	155	0
	Recheck	2	1
	Open	1	33


Total of 192 Cows Tested

Sensitivity* = 99.4% (95% CI: 96.0-100%)

Specificity* = 100% (95% CI: 87.3-100%)

**Recheck results excluded*

Assay Response of Individual Cows



IDEXX Milk Pregnancy Test

Results – Herd B

Midwestern Dairy: 120 cows sampled

All cows > 60 DIM and not bred or >40 days post insemination

Ultrasound or Palpation

**IDEXX
Milk ELISA**

	Pregnant	Open
Pregnant	104	0
Recheck	0	1
Open	3	12


Total of 120 Cows Tested

Sensitivity* = 97.2% (95% CI: 91.6-99.4%)

Specificity* = 100% (95% CI: 71.3-100%)

**Recheck results excluded*


Assay Response of Individual Cows



IDEXX Milk Pregnancy Test

Detection of Open Cows

- 290 cows from 5 geographically diverse herds
- All cows > 60 DIM and not bred or >40 days post insemination
- Bred cows confirmed open by ultrasound or palpation



- Overall specificity = 96.7% (95% CI: 94.9-98.9%)
- Fewer open samples available in late gestation, but specificity is very good (100% after day 90)
- Specificity is lower prior to day 40 of pregnancy, likely due to high rates of embryonic death and residual PAGs in the cow's system

Integrating Milk Pregnancy Testing into Routine Herd Recording Programmes

- Field
 - Versatile Application
 - ▶ Early detection and confirmations
 - ▶ Preselected by days since bred

- Laboratory
 - Simple ELISA platform
 - Uses existing samples

- Data Processing Centers
 - Rapid turnaround
 - Reproduction reports
 - Hot sheet



025 - VET - COW PREG CK, POST PARTUM, OPEN, NOT BRED, & FRESH

UNIV OF IL DAIRY HERD - 33090003
Ref:06-17

G	Reason	Days	Svc	C	Bred	Days	Most	Recent	Health	Code	BCS	
R	Barn	on	In	Heat	No	D	Heat	Sinc	Current	Lactation	Mst	
P	Name	Report	Milk	Int	Br	E	Date	Br-H	MMDDYY	Code T	Remarks	Rct
3	DARLENE	PREGCK	332	31	7		5-09	40	041304	HRMX Y	POSTSYNCH	GNRH
3	LASSY	PREGCK	168	22	3		5-07	42	042004	HRMX Y	POSTSYNCH	GNRH
3	LEOTA	PREGCK	358	42	7	R	5-06	43	060804	PALP N	PRH	EED?
3	EDNA	RECK-P	203	56	2	P	3-25	85	042704	PALP N	PRH	
3	MILLY	PREGCK	297	42	6	R	5-06	43	060804	HRMX Y	OVSYNCH	GNRH1
3	BEAUTY	PREGCK	430	49	7	R	5-06	43	060804	PALP N	PRH	
3	CANDEE	PREGCK	108		1	R	5-06	43	060804	PALP N	PRH	
3	DARLENE	RECK-P	269	42	4	P	4-01	78	050404	PALP N	PRH	3.5
3	7453	RECK-P	188	42	2	P	3-25	85	042704	PALP N	PLH	
3	SUZANN	RECK-P	155		1	P	3-18	92	042704	PALP N	PRH	
3	CHEERFU	PREGCK	187	42	3	R	5-06	43	060804	PALP N	PRH	
3	PRECISI	RECK-P	204	49	2	P	3-18	92	042004	PALP N	PLH	
3	SUZANNA	RECK-P	145		1	P	3-25	85	050404	PALP N	PLH	
3	LASSIE	PREGCK	201	20	4		5-09	40	041304	HRMX Y	OVSYNCH	GNRH1

Conclusions

- PAGs can be detected by the IDEXX Milk Pregnancy Test throughout gestation
- The IDEXX Milk Pregnancy Test offers a expanded testing options for herd recording laboratories
- Herd recording samples for pregnancy diagnosis - hassle-free and cost-effect reproductive management in dairy herds
- Application of a milk-based pregnancy test coupled with management data from herd recording could help to improve reproductive efficiency in dairy herds



Test With Confidence™



For more information please contact...

Hannah Pearse

Marketing Manager, IDEXX Inc.

Email: hannah-pearse@idexx.com

Phone: 0044 (0) 7584 681 031