




# Innovative analytical and automation solutions to extract the maximum value of every milk samples for optimal herds management

By Pierre Broutin  
Managing, R&D Director/Western Europe  
[pbroutin@bentleyinstruments.eu](mailto:pbroutin@bentleyinstruments.eu)  
[www.bentleyinstruments.eu](http://www.bentleyinstruments.eu)



# Introduction

---

EU milk production is shifting to a more sustainable model associated with the reduction of milk production over the coming years.

This shift will probably keep increasing market segmentation with consumers demand for high-quality, healthy and sustainable dairy products, especially for organic, low-fat, low-sugar, no allergens, probiotic and fortified dairy products. Cheese and whey are expected to benefit the most from this shift.

Sustainability is increasingly important to the dairy industry. Producers will need to adapt to reduce their carbon footprint and environmental impact, while maintaining their business profitability.

How to produce high quality dairy products while taking into account all new societal and environmental demands in terms of health, sustainability, animal welfare, carbon footprint, at a price accepted by consumers? What are the new milk components and functionalities of interest? Can milk composition and functionalities be optimized to meet markets new demand?


Reference methods are typically time consuming and quite expensive. New fast alternative, standardized and cost effective analytical methods are necessary to promote and monitor these new compositional and functional evolutions throughout the dairy chain.

FT-MIR technology is one of the obvious candidate being already used since 1990 on a large scale in the dairy industry and central milk testing laboratories for milk payment, dairy herd management and incoming milk and process monitoring.


# Innovative Analytical & Automation Solutions for the Dairy Industry

## Bentley Instruments, Inc (USA)


- Company established in 1983 (near Minneapolis)
- Presents in > 60 countries
- > 110 employees
- 5 Subsidiaries in Europe + distribution network
- All our instruments IDF/ISO/ICAR certified
- Market: dairy industry exclusively



# Solutions for the Dairy Plants for the rapid and highly accurate determination of milk chemical and hygienic composition



**DairySpec FT**  
Milk/Dairy Products  
Up to 64 chemical components  
ICAR certified



**SomaCount FC**  
Somatic Cells  
ICAR certified




**New BactoCount IBCM3**  
Total bacteria/Somatic cells  
ISO 16140 certified  
**By far, the fastest method on the market for total bacteria testing**






# Solutions for the Central Milk Testing Laboratories

## BactoCount: Rapid and highly accurate determination of milk hygienic composition



**New BactoCount IBC 200 Multiplex**  
**Total Bacteria/Somatic Cells**  
**Up to 200 samples/hour**  
**ISO 16140 Microval Certified**  
**Cow/Sheep/Goat**



**MICROVAL<sup>®</sup>**   
**NEI**

**CERTIFICATE OF COMPLIANCE**  
**LLOYD'S REGISTER QUALITY ASSURANCE**

hereby declares that the certification assessment has demonstrated that

**BactoCount IBC & IBCm**

Manufactured and supplied by:  
 Bentley Instruments  
 4004 Peavey Road, Chaska  
 MN 55138 USA

has been validated and revealed to be at least equivalent to the reference methods as demonstrated by the validation study report. The summary of the validation report is available on the MicroVal website: [www.microval.org](http://www.microval.org)

- Reference methods:
1. EN-ISO 4833-1:2013 Microbiology of the food chain – Horizontal method for the enumeration of microorganisms – Part 1: Colony count at 30 degrees C by the pour plate technique.
  2. EN-ISO 4833-2:2013 Microbiology of the food chain – Horizontal method for the enumeration of microorganisms – Part 2: Colony count at 30 degrees C by the surface plating technique.

Scope: raw cow's milk

The validation and certification has been performed in accordance with EN ISO 16140-2:2016, the EURL MMP document - Validation criteria of instrumental methods for enumeration of total flora in raw milk, version 2, 21 December 2011 and the MicroVal Rules and Certification Scheme version 8.

Certificate no.: 2013LR44

First approval date: 16 December 2014

Renewal date: 03 March 2019

Expiry date: 15 December 2022

ISSUED BY: Lloyd's Register Nederland B.V.  
 Rotterdam, The Netherlands

Certificate no.: 2013LR44


11-03-2019

Page 1 of 1

K.P. v.d. Mandelelaan 41a, 3062 MB Rotterdam, The Netherlands. KvK nr.: 24247948  
 This approval is carried out in accordance with the LQA assessment and certification procedures and monitored by LQA.



**BactoCount IBC 50-150**  
**Total Bacteria**  
**ISO 16140 Microval Certified**



# Solutions for the Central Milk Testing Laboratories

Only Combi systems ICAR certified for cow, goat and sheep milk testing



**CombiFTS 400/500/600**  
FTIR+FC (up to 64 components)



**IR spectra standardized in real time  
without any reagent**



**DairySpec Combi 100/200/300**  
FTIR+FC (up to 64 components)

## Bentley Combi FTS 600 & IlaS 3000 robot

Fully automated and standardized samples preparation for the highly accurate determination of **milk chemical and somatic cells composition**



### Main characteristics :

- Up to 600 samples/hour (Combi FTS)
- Up to 64 chemical parameters analyzed simultaneously
- Spectra standardized\* in real time (without reagents)
- Samples handling fully automated and standardized for optimal analytical performances

### The samples are automatically:

- Identified (RFID, code-barres...)
- Heated to 40°C or Cooled at 4°C
- Inverted and uncapped
- Analyzed
- Replaced in their original position

\*patented





# Quantitative IR analysis

Fatty Acids origin profiling to feed dairy herds for optimal milk fat composition

**Milk contains more than 400 individual fatty acids:**  
(over 100 publications over last decade)

## De Novo/Short Chain Fatty Acids (18-30%)

C4:0, C6:0, C8:0, C10:0, C12:0, C14:0, C14:1

## Preformed/Long Chain Fatty Acids (30-45%)

C18:0, C18:1, C18:2

## Mixed/Medium Chain Fatty Acids

C16:0, C16:1

## Traditional Fatty Acids

Saturated, Unsaturated, Mono & Polyunsaturated, Oleic (C18:1),  
Palmitic (C16:0) and Stearic (C18:0)



MOATE ET AL.

**Table 2.** The proportions of 26 individual fatty acids (mg/g) in the total milk fatty acids from milks described in 28 publications

Fatty acid <sup>1</sup>	n <sup>2</sup>	Mean ± SD	Median	Minimum	Maximum
4:0	95	31.3 ± 6.8	29.9	18.4	49.2
6:0	111	19.4 ± 5.2	19.2	6.3	32.3
8:0	111	11.7 ± 3.5	12.0	4.8	20.9
10:0	111	24.8 ± 7.3	25.0	10.3	39.4
12:0	111	29.9 ± 8.5	29.0	15.0	52.2
14:0	111	103.8 ± 17.1	101.7	63.3	135.0
14:1c9	101	10.8 ± 3.6	10.9	3.5	21.3
15:0	88	10.5 ± 3.3	10.2	4	22.6
16:0	120	285.1 ± 49.8	281.5	147.1	462.1
16:1c9	109	17.3 ± 6.3	17.3	4	36.5
C <sub>17</sub>	78	7.3 ± 3.5	5.9	3.3	16.7
18:0	120	105.1 ± 35.9	99.7	30.6	268.7
18:1t6-8	33	4.6 ± 2.1	4.9	1.2	9.6
18:1t9	37	4.4 ± 2.0	4.4	1.4	11.4
18:1t10	30	13.1 ± 15.2	8.1	0.3	64.7
18:1t11	90	33.3 ± 21.8	32.6	5.8	99.5
18:1t12	19	6.5 ± 3.6	6.3	0.9	12.7
18:1c9	120	205.0 ± 53.5	199.6	70.3	371.4
18:2c9,c12	120	31.3 ± 21.1	26.6	5.1	133.0
18:2c9,t11	76	10.2 ± 6.0	8.4	2.8	24.5
18:2t10,c12	35	0.4 ± 0.3	.3	0	1.4
18:2c11,t13	6	0.4 ± 0.3	0.3	0.2	0.9
OCLA	25	1.5 ± 1.4	1.1	0	4.2
18:3	114	5.9 ± 3.6	4.9	0.2	19.0
20:0	30	1.5 ± 0.6	1.5	0.4	3.0
20:5	39	1.0 ± 1.1	0.5	0	4.8
22:6	31	0.7 ± 0.7	0.4	0	2.6
Others	120	75.1 ± 56.2	76.4	0	237.3
Total CLA	82	10.3 ± 6.6	8.3	3	28.4
Total 18:1 trans	94	42.5 ± 26.3	39.5	8.6	145.1
Total de novo	120	232.6 ± 42.4	237.9	136.6	300.6
Total C <sub>16</sub>	120	300.9 ± 52.7	303.0	154.2	462.1
Total preformed	120	466.5 ± 75.8	477.4	315.6	641.3


<sup>1</sup>c = cis; t = trans; OCLA = other conjugated linoleic acids; CLA = conjugated linoleic acid; total 18:1 trans = sum of 18:1t6-8, 18:1t9, 18:1t10, 18:1t11, and 18:1t12 isomers; total de novo = sum of 4:0 to 15:0 fatty acids; total C<sub>16</sub> = sum of 16:0 and 16:1; total preformed = sum of all milk fatty acids with more than 17 carbon atoms.

<sup>2</sup>n = total number of dietary treatments contributing to each mean.





# New biomarkers for ketosis detection Blood BHB\* (from milk spectrum)



(19) (11) EP 3 058 341 B1

(12) FASCICULE DE BREVET EUROPEEN

(45) Date de publication et mention de la délivrance du brevet: 16.05.2018 Bulletin 2018/20

(21) Numéro de dépôt: 14824867.7

(22) Date de dépôt: 17.10.2014

(51) Int. Cl.:  
G01N 21/27 (2006.01) G01N 21/35 (2014.01)  
G01N 33/04 (2006.01) G01N 33/49 (2006.01)  
A61B 5/00 (2006.01) G01N 21/31 (2006.01)  
C12Q 1/32 (2006.01) G01N 33/50 (2006.01)  
G01N 21/3577 (2014.01) G01N 33/487 (2006.01)

(86) Numéro de dépôt international: PCT/FR2014/052650

(87) Numéro de publication internationale: WO 2015/055966 (23.04.2015 Gazette 2015/16)

(54) DETERMINATION DE LA CONCENTRATION D'UN COMPOSANT DANS UN FLUIDE D'UN ANIMAL PAR ANALYSE SPECTROSCOPIQUE D'UN AUTRE FLUIDE  
BESTIMMUNG DER KONZENTRATION EINER KOMPONENTE IN EINEM FLUID EINES TIERES DURCH SPEKTRALANALYSE EINES ANDEREN FLUIDS  
DETERMINATION OF THE CONCENTRATION OF A COMPONENT IN ONE FLUID OF AN ANIMAL BY SPECTROSCOPIC ANALYSIS OF ANOTHER FLUID

## Prediction of blood components from milk spectra based on metabolic disorders modeling

59000 Lille (FR)

(74) Mandataire: Argyma  
36, rue d'Alsace Lorraine  
31000 Toulouse (FR)

(56) Documents cités:  
• HANSEN ET AL: "Screening of Dairy Cows for Ketosis by Use of Infrared Spectroscopy and Multivariate Calibration", JOURNAL OF DAIRY SCIENCE, AMERICAN DAIRY SCIENCE ASSOCIATION, US, vol. 82, no. 9, 1999, pages 2005-2010, XP026993664, ISSN: 0022-0302 [extrait le 1999-09-01]

[extrait le 2007-04-01]  
• OTO HANUS ET AL: "Milk acetone determination by the photometrical method after microdiffusion and via FT infra-red spectroscopy", JOURNAL OF AGROBIOLOGY, vol. 28, no. 1, 2011, XP055119188, ISSN: 1803-4403, DOI: 10.2478/v10146-011-0004-9

• ROBERT P ET AL: "MULTIVARIATE ANALYSIS APPLIED TO NEAR-INFRARED SPECTRA OF MILK", ANALYTICAL CHEMISTRY, AMERICAN CHEMICAL SOCIETY, US, vol. 59, no. 17, 1987, pages 2187-2191, XP000857119, ISSN: 0003-2700, DOI: 10.1021/AC00144A038


Il est rappelé que: Dans un délai de neuf mois à compter de la publication de la mention de la délivrance du brevet européen au Bulletin européen des brevets, toute personne peut faire opposition à ce brevet auprès de l'Office européen des brevets, conformément au règlement d'exécution. L'opposition n'est réputée formée qu'après le paiement de la taxe d'opposition. (Art. 99(1) Convention sur le brevet européen).

EP 3 058 341 B1

# Milk Amyloid A (MAA): a potent biomarker for early mastitis diagnosis and milk quality monitoring

Natural Clinical Mastitis		
Sensitivity	Specificity	Reference
93.0%	100%	Eckersall et al. 2001

Natural Subclinical Mastitis		
Sensitivity	Specificity	Reference
90.6%	98.3%	Safi et al. 2009
92.3%	92.1%	Shirazi-Beheshtiha et al. 2012



MAA is the only Acute Phase Protein produced directly by the epithelium of the udder in response to bacterial infection in the udder and as such is an immediate and direct marker of infection



Exclusive worldwide distributor



# Bentley Combi FTS

## New potential developments: Milk-based Genetic analysis & selection







### Sequence-based genome-wide association study of milk mid-infrared wavenumbers in dairy cattle

Kathryn Tiplady

T. Lopdell<sup>1</sup>, E. Reynolds<sup>2</sup>, R. Sherlock<sup>1</sup>, M. Keehan<sup>1</sup>, T. Johnson<sup>1</sup>, J. Pryce<sup>3,4</sup>, H. Blair<sup>2</sup>, S. Davis<sup>1</sup>, M. Littlejohn<sup>1,2</sup>, D. Garrick<sup>2</sup>, R. Spelman<sup>1</sup> and B. Harris<sup>1</sup>

<sup>1</sup>Livestock Improvement Corporation, Hamilton, NZ  
<sup>2</sup>Massey University, AL Rae Centre, Hamilton, NZ  
<sup>3</sup>Agriculture Victoria, AgriBio, VIC, Australia  
<sup>4</sup>La Trobe University, VIC, Australia



Bonferroni-threshold:  $0.01/17.7m/895=6.3 \times 10^{-13}$ .


Significant associations conserved across multiple wavenumbers: DGAT1 significant for over 750 wavenumbers!



## Development of new IR milk phenotypes necessary to keep improving dairy herds management and sustainability

- **Feeding optimization:** fatty acids profile, urea, protein...
- **Early mastitis detection:** SCC, Milk Amyloid A, lactoferrin, citrate ...
- **Early Ketosis detection:** Blood BHB\*, ketone bodies, citrate, fatty acids profile ...
- **Breeding/Genetic selection:** Standardized IR spectra\*, fat & protein profiling, urea, ketone bodies...
- **Environmental impact mitigation:** nitrogen, methane emission (from fatty acids profile), phosphorus....



\*patented




# A new holistic approach for the optimum valorization of every milk samples

## Why?

- **Streamline mass & multiplex testing at the sample level** to reduce operation costs while extracting the maximum value of every samples of interest by intelligently combining multiple analytical methods
- **Standardize milk samples preparation and handling before testing**, a critical step to secure analytical methods standardization and performances for worldwide results equivalence
- **Increase laboratory throughput** with potential speed up to 600 samples/hour and deal more efficiently with peaks in samples distribution with samples automatically dispatched to open analyzers



# Innovative analytical and automation solutions to extract the maximum value of every milk samples for optimal herds management







## Conclusion

---

The development of a holistic model combining latest high throughput standardized methods in combination with complete laboratory automation makes it possible to streamline mass & multiplex testing at the sample level for the benefit of the dairy industry.



---

Thank you for your attention!  
Looking forward to meeting you on our stand

[pbroutin@bentleyinstruments.eu](mailto:pbroutin@bentleyinstruments.eu)

[www.bentleyinstruments.eu](http://www.bentleyinstruments.eu)

