

Use of a highly accurate enzymatic method to evaluate the relationship between Milk Urea Nitrogen, milk composition and yield on bulk and individual milk samples

Rapid determination of somatic cells and total flora in cow, goat, sheep and buffalo milk

By Pierre Broutin
Senior Scientist/Application Development
Managing Director/ Western Europe
Bentley Instruments
pbroutin@bentleyinstruments.com

THE AS62 DHIA LABORATORY UREA EXPERIENCE

MATERIALS AND METHOD

- > Started a UREA testing program in 1995 with infrared technology. Accuracy was found to be insufficient to explain the observations made on the field and fully optimize the feeding program
- > Decided to switch to chemical analysis in 1998 with the introduction of the ChemSpec enzymatic method
- > Retrospective Analysis of individual cows records (2003-2005) from monthly DHIA tests collected by the AS62 DHIA laboratory (Marceuil, France)
- > ~ 600 000 individual cow milk samples tested for fat, protein, somatic cells and urea
- > Determination of the fat, protein and somatic cells content on the Bentley B2000 Milk Analyzer
- > Determination of the Somatic Cells count on the Bentley Somacount (flow cytometry)
- > Determination of the Urea Content on the Bentley ChemSpec (enzymatic)
- > HERD Mean Values grouped by Urea Classes (1 mg/dl increment)
- > Develop their own specific diagnostic grid to optimize the use of urea values

MATERIALS AND METHOD

ChemSpec 150 Milk Urea Analyzer



C. LHOMME – AS62 General Manager
G. HEUMEZ – AS62 Technical Director

- 600 000 milk samples tested for urea every year
- 4 ChemSpec 150

ChemSpec Method Principle

- > Enzymatic reaction + colorimetric detection (green)
- > Differential measurement (with/without enzyme)
- > 150 samples/hour
- > Fully automated
- > No sample preparation (4°C - 40°C)
- > Method highly specific → One Calibration Std.
- > Method highly accurate ($s_{y,x} < 1$ mg/dl urea)

ChemSpec 150 Accuracy
Bulk Milk Samples

Instrument	N	Sy,x	Reference Methods	Sources
ChemSpec 150	10	0.17	Enzymatic + pH Differential	Cecalait Ring test (1997)
ChemSpec 150	10	0.19	AFNOR NF V04-217	Cecalait Ring test (1997)
ChemSpec 150	10	0.21	DMAB	Cecalait Ring test (1997)
ChemSpec 150	10	0.26	Enzymatic + Colorimetric	Cecalait Ring test (1997)
ChemSpec 150	78	0.77	DMAB	Broutin(1998)
ChemSpec 150	30	0.88	Enzymatic + pH Differential	Broutin(1997)
ChemSpec 150	233	0.96	Enzymatic + pH Differential	Associazione Regionale Allevatori Laboratory, Crema, Italy (2000)
ChemSpec 150	96	0.74	Enzymatic + pH Differential	Compilation of the National DHIA MUN reports (Jan to August 2002)
ChemSpec 150	49	1.58	Enzymatic AFNOR NF V04-217	Cecalait Evaluation (1999)
ChemSpec 150		<1.00	Enzymatic + pH Differential	Manufacturer specification
Filter IR		<1.5		Manufacturer specifications
FTIR		<1.1		Manufacturer specifications

ChemSpec Accuracy >> FTIR and Filter IR Methods


ChemSpec 150 Accuracy
Individual Milk Samples

Instrument	N	Sy,x	Reference Methods	Sources
ChemSpec 150	129	0.82	Enzymatic + Conductivity	Milchprüfung Niederosterreich Laboratory, Austria (2000)
ChemSpec 150	139	1.00	DMAB	Broutin(1998)
ChemSpec 150	158	1.26	Enzymatic + pH Differential	Broutin(1999)
ChemSpec 150	98	2.10	Enzymatic AFNOR NF V04-217	Cecalait Evaluation (1999)
ChemSpec 150		<1.5	Enzymatic + pH Differential	Manufacturer specification
Filter IR		<1.5		Manufacturer specifications
FTIR		<1.1		Manufacturer specifications

ChemSpec Accuracy >> FTIR and Filter IR Methods

THE AS62 DHIA LABORATORY UREA EXPERIENCE

HERD UREA FREQUENCY DISTRIBUTION




~ 23000 Herds ~ 600 000 samples tested over 3 years

35th ICAR Session KUOPIO, FINLAND 2006

MUN = UREA/2.145

Relationship between Milk Urea and Fat Content (~23000 herds)




UREA CLASSES (MG/DL)	MEAN FAT (GL)
<10	41,75
>10 <20	41,09
>20 <30	40,31
>30 <40	39,94
>40 <50	39,83
MAX DIFF	1,92

→ STRONG NEGATIVE CORRELATION BETWEEN UREA AND FAT CONTENT

35th ICAR Session KUOPIO, FINLAND 2006

Relationship between Milk Protein and Urea content (~23000 herds)




UREA CLASSES (MG/DL)	MEAN TP (G/L)
<10	31,25
>10 <20	31,81
>20 <30	31,94
>30 <40	32,00
>40 <50	32,06
MAX DIFF	-0,82

Source: AS62

→ Positive correlation between urea and protein content

35th ICAR Session KUOPIO, FINLAND 2006

Relationship between Milk Somatic Cells and Urea content (~23000 herds)




Source: AS62

→ negative correlation in lower urea range (<30 mg/dl)
→ positive correlation in higher urea range (> 30 mg/dl)

35th ICAR Session KUOPIO, FINLAND 2006

Relationship between Milk Production Level and Urea content (~23000 herds)




Source: AS62

→ Strong positive correlation (< 30 mg/dl)
→ Negative correlation (> 40 mg/dl)

35th ICAR Session KUOPIO, FINLAND 2006

UREA CLASSES	MILK YIELD
<10	7399
>10 <20	8323
>20 <30	8885
>30 <40	9030
>40 <50	8803
MAX DIFF	1631

Effect of the Milk Urea Monitoring on Milk Production Level




Source: AS62

MONTH	WITH	WITHOUT	DIFF	%
01	8287	8619	-331,87	-4,0
02	8518	8929	-411,09	-4,8
03	8879	9354	-474,72	-5,3
04	9520	9922	-402,01	-4,2
05	9574	9957	-383,21	-4,0
06	9160	9512	-352,66	-3,9
07	8990	9325	-335,35	-3,7
08	8908	9248	-340,23	-3,8
09	8448	8741	-293,63	-3,5
10	7976	8347	-370,73	-4,6
11	8124	8373	-249,24	-3,1
12	8252	8507	-255,34	-3,1
MEAN	8689	9053	-364	-4,19

35th ICAR Session KUOPIO, FINLAND 2006

Effect of Milk Urea Monitoring on Milk Protein Content




Source: AS62

MONTH	WITH	WITHOUT	DIFF	%
01	31,79	31,85	-0,06	-0,19
02	31,66	31,70	-0,04	-0,12
03	31,55	31,67	-0,12	-1,03
04	31,64	32,06	-0,42	-1,34
05	31,69	32,13	-0,44	-1,39
06	31,33	31,72	-0,38	-1,23
07	31,66	31,74	-0,18	-0,58
08	31,30	32,03	-0,73	-2,16
09	31,68	32,03	-0,35	-1,11
10	31,39	31,62	-0,23	-0,73
11	32,20	32,25	-0,05	-0,17
12	32,10	32,25	-0,15	-0,47
MEAN	31,68	31,97	-0,29	-0,92

35th ICAR Session KUOPIO, FINLAND 2006

Effect of Urea Monitoring on Somatic Cells Level



MONTH	WITH	WITHOUT	DIFF	%
01	329	308	21,00	6,79
02	313	303	10,00	3,30
03	306	310	-4,24	-1,38
04	342	312	29,33	9,41
05	320	295	24,98	8,47
06	331	314	16,85	5,37
07	313	295	16,75	5,68
08	310	288	21,58	7,49
09	273	262	10,84	4,14
10	261	250	11,08	4,43
11	284	282	1,87	0,66
12	305	302	2,74	0,91
MEAN	307	290	17	5,70

35th ICAR Session KUOPIO, FINLAND 2006

ROCKY MOUNTAIN DHIA EXPERIENCE


MATERIALS AND METHOD

- 24 months retrospective analysis of individual cows records from monthly DHIA tests collected by Rocky Mountain DHIA
- 49 Holstein Herds
- 69724 individual cow records (fat, protein, somatic cells and urea)
- Determination of the fat, protein and somatic cells content on the *Bentley B2000* Milk Analyzer
- Determination of the Somatic Cells count on the *Bentley Somacount* (flow cytometry)
- Determination of the Urea Content on the *Bentley ChemSpec*
- **INDIVIDUAL COW VALUES** grouped by Urea Classes (~4 mg/dl increment)

35th ICAR Session KUOPIO, FINLAND 2006


The Association between milk urea nitrogen and DHI production variables in Western Commercial Dairy Herds

(Urea values data treated at the **COW** Level)



J. Dairy Sci. 86:3008-3015 by R.G. Johnson, A.J. Young

The Association between milk urea nitrogen and DHI production variables in Western Commercial Dairy Herds



35th ICAR Session KUOPIO, FINLAND 2006


CONCLUSIONS

THESE STUDIES SUGGEST THAT:

- THERE ARE STRONG RELATIONSHIPS BETWEEN UREA, FAT, PROTEIN AND SOMATIC CELLS CONTENT AS WELL AS MILK PRODUCTION.
 - THE ACCURATE DETERMINATION OF UREA ON INDIVIDUAL COWS CAN BE A VERY USEFUL TOOL TO OPTIMIZE FEEDING PROGRAMS AND INCREASE MILK PRODUCTION AND PROTEIN CONTENT
- A VERY SUCCESSFUL EXPERIENCE:**
- THE AS62 DHIA LABORATORY IS NOW PERFORMING 600 000 UREA TESTS ON INDIVIDUAL COWS EVERY YEAR (60% OF THE PRODUCERS) ON THE CHEMSPEC (4)
 - THE NUMBER OF TESTS IS STILL INCREASING DUE TO THE TEST ADDED VALUE (ONLY COMPONENT)
 - THE MARKET DEMAND WILL ORIENTATE OUR CHOICES BUT THE CHEMSPEC SPEED COULD BE SIGNIFICANTLY INCREASED TO REACH OUR CUSTOMERS NEEDS

35th ICAR Session KUOPIO, FINLAND 2006

Rapid Determination of Somatic Cells in Goat Milk by Flow Cytometry CECALAIT Evaluation



Source: CECALAIT LETTER N°54

Classes	0-2000	0-1000	1000-2000
N	44	24	20
Sr	13	9	17
Sr(%)	1,40%	1,30%	1,30%

Source: CECALAIT LETTER N°54

- Somacount Somatic Cells Counter (FC) calibrated with CECALAIT COW CRM samples
- Sy,x (%) = 4,7%
- Sr(%) = 1,4%

35th ICAR Session KUOPIO, FINLAND 2006

Rapid and accurate determination of milk total flora by flow cytometry BactoCount IBC & IBCm




BactoCount 50-150



BactoCount IBCm

35th ICAR Session KUOPIO, FINLAND 2006

BactoCount IBC Universal Calibration Cow, Sheep, Goat and Buffalo Milk



35th ICAR Session KUOPIO, FINLAND 2006

THANK YOU FOR YOUR ATTENTION

pbrouin@bentleyinstruments.com

35th ICAR Session KUOPIO, FINLAND 2006