

# Specification of straw color and identification

Description of Swissgenetics' specification of straw color and identification, effective since 16 March 2011.

# 1 Straw specification

# 1.1 Size

French straws 0.25 ml

# 1.2 Colors

Swissgenetics has defined a specific straw color for each principal breed in Switzerland. Additional colors are used for special semen straws (articles).

Breed	Code	Semen straw color	
Brown Swiss	BS	green (transparent)	
Swiss Fleckvieh SF		red (transparent)	
Holstein	HO	grey (transparent)	
Limousin	LM	ocher	
Original Braunvieh	OB	green (transparent)	
Red Holstein	RH	red (transparent)	
Simmental	SI	red (transparent)	

The color for other (exotic) breeds is not defined and is chosen individually for each semen production.

۰.					
	Product	Semen straw color			
	Heterospermic doses	orange			
	Sexed semen (X)	lilac (Cryovet 662) (transparent)			
	Sexed semen (Y)	blue (Cryovet 262) (transparent)			
	others	light			

The color for special articles is defined as follows:

# 2 Straw identification

# 2.1 Sequence of the imprinting on the straw

	Production date		Article code	Batch number		Name of the bull	Breed code
CHAGAI02B	140514		N	046	CH 120.0917.3984.8	HARRY	BS
19608425046							

## 2.2 Detailed information

#### EU code

Since January 1st, 2007, the official code for the semen collection center Mülligen assigned by the Swiss Federal Food Safety and Veterinary Office is CH-AG-AI-02B (before Code: CH-AI-2B). In order to save



Page 2/3

space, the EU code is printed without hyphens.

# **Production date**

The production date is printed as DDMMYY<sup>1</sup>, with

- DD =day
- MM = month
- YY = last two digits of the year of production (e.g. 14 for 2014)

# Barcode (128C)

The barcode is a machine-readable information according to the ICAR Guidelines. It is printed in type 128C format and consists of the following elements:

NAAB stud code	196: Registration number of the semen collection center Mülligen at NAAB
Internal number Bull	Unique number of Bull x Article (5 digits), as assigned by the internal central database
x Article	of Swissgenetics
Batch number	Number within Bull x Article as defined below (3 digits)

Examples: see Appendix

Please note: A barcode of type 128C includes an internal checksum digit.

#### Article code

This code represents the specificity of the semen doses, e.g. sexed semen. Different products (articles) can be produced from the same bull, depending on the processing procedure. In order to be able to distinguish these products, a process-specific code is necessary.

Code	Specification
Ν	Standard/conventional semen
Н	Heterospermic doses
S	SpermVital
Х	Sexed semen (X)
Y	Sexed semen (Y)

#### Batch number

A batch is in the terminology of Swissgenetics a certain number of semen straws of an article processed the same way, which origin from one or more ejaculates of a bull at the same day. Batch numbers are assigned to ejaculates processed **within Bull x Article**. If at the same day two ejaculates of a bull are merged for producing one kind of article, one batch number is assigned to the mixture. If the ejaculates are processed separately, different batch numbers are assigned.

#### Animal ID

Official herdbook identification of the bull

#### Name of the bull

Name (short form) of the bull (maximum of 12 characters)

#### **Breed code**

Breed code according to herdbook registration

<sup>&</sup>lt;sup>1</sup> From March 16<sup>th</sup>, 2011 till March 31<sup>st</sup>, 2014, the production date was printed as YYDDD with DDD = n<sup>th</sup> day within year (e.g. 054 for February 23<sup>rd</sup>)



# Page 3/3

# Appendix: Barcode example of a single bull producing different articles

Production date	Article	Barcode Type 128C	In clear text
09 Jan 2017	Conventional		<mark>196</mark> 00430 <mark>153</mark>
12 Jan 2017	Conventional		<mark>196</mark> 00430 <mark>154</mark>
16 Jan 2017	Conventional		<mark>196</mark> 00430 <mark>155</mark>
19 Jan 2017	Conventional		<mark>196</mark> 00430 <mark>156</mark>
24 Jan 2017	SpermVital		<mark>196</mark> 11824 <mark>013</mark>
26 Jan 2017	Sexed (X)		<mark>196</mark> 07966 <mark>095</mark>
26 Jan 2017	Sexed (X)		<mark>196</mark> 07966 <mark>096</mark>
02 Feb 2017	Sexed (X)		<mark>196</mark> 07966 <mark>097</mark>
02 Feb 2017	Sexed (X)		<mark>196</mark> 07966 <mark>098</mark>
06 Feb 2017	Conventional		<mark>196</mark> 00430 <mark>157</mark>
09 Feb 2017	Conventional		<mark>196</mark> 00430 <mark>158</mark>
13 Feb 2017	Conventional		<mark>196</mark> 00430 <mark>159</mark>
20 Feb 2017	Conventional		<mark>196</mark> 00430 <mark>160</mark>
23 Feb 2017	Conventional		<mark>196</mark> 00430 <mark>161</mark>
27 Feb 2017	Conventional		<mark>196</mark> 00430 <mark>162</mark>
06 Mar 2017	Conventional		<mark>196</mark> 00430 <mark>163</mark>
13 Mar 2017	Conventional		<mark>196</mark> 00430 <mark>164</mark>
15 Mar 2017	Conventional		<mark>196</mark> 00430 <mark>165</mark>
23 Mar 2017	Sexed (X)		<mark>196</mark> 07966 <mark>099</mark>
28 Mar 2017	Conventional		<mark>196</mark> 00430 <mark>166</mark>
03 Apr 2017	Conventional		<mark>196</mark> 00430 <mark>167</mark>
05 Apr 2017	SpermVital		<mark>196</mark> 11824 <mark>014</mark>

## where

<mark>196</mark> 00430/11824/07966 <mark>153</mark> ... = NAAB stud code of the semen collecting center Mülligen

= unique internal number of Bull x Article

= batch number within Bull x Article