

Perspectives for proactive dairy herd management based on on-farm milk analysis- and potential implications for future milk recording

By M. Sc. Tove Asmussen and M. Sc. Ph D Carsten Ridder, Lattec

lattec

Purpose of this presentation

- To present some results from the Biosens models, on DCRC data and from the two commercial test farms
- To discuss the need for ICAR specifications to on-farm measurements


lattec

lattec

- Established in 2001
- 50:50 partnership between DeLaval and FOSS
- Development, market preparation and production
- Sales via owners:
- DeLaval to the farmers
- FOSS via existing customers
- Situated in Hillerød

lattec

Lattec's solution, an overview...



lattec

Standard Operating Procedures are created!


- Biosens, a research consortium between DIAS, DCRC and Lattec
- Biosens builds models, Lattec implements, - and we test the models together, creating SOP's.

We focus at:

- Mastitis
- Heat, reproductive disorders and pregnancy
- Ketosis
- Imbalances in feeding


lattec

Mastitis, - typical example




lattec

Mastitis: Alert Nov 9th, detected Nov 12th, treated Nov 15th (com. test farm)




lattec

Ketosis, Lattec alert Oct 19th and Nov 1st, owner detection Oct 24th and Nov 4th



lattec

Reproduction



lattec

Energy model, - including fat and protein?

- Lattec's purpose for measuring fat and protein
 - To detect upcoming metabolic disorders
 - To improve feeding and reproduction management
- ICAR's purposes for measuring fat and protein
 - For improving feeding and reproduction management
 - For breeding purposes

So we do have mutual interests!

lattec

Mutual target!

- Holistic view. Not only talk repeatability and bias
- Go for an overall improvement of the whole system. Not only parts of it.


lattec

Contributors to inaccuracy in measured results

- Correct sample ID
- Sample quality at measurement
- Measurement frequency
- Repeatability and accuracy (bias) of the measurement.
- Day-to-day variation in fat/protein content of the milk sample

lattec

Day to day variation in fat content (DK results)




lattec

Day-to-day variation, -now two data sets:

- **DK-dataset:**
 - 1 herd, 11 cows randomly selected, monitored all lactation
- **CA-dataset:**
 - 14 herds, 10 cows from each, each cow is monitored for 5 days.

lattec

Standard deviation (fat) over 5 days



Average CV (DK): 5.25


Average CV (CND): 6.40

Average CV Parlor 1: 6.43

Average CV Parlor 2: 6.54


lattec

DHIA today



lattec

How do these results compare to today's practice?



lattec

One approach....

$$CV\% = (SD \times 100\%) / \text{av fat } \%, - \text{ Av Fat}\% = 4.5\%$$

$$CV (MSC) = 0.5 \%, \text{ or } SD(\text{instrument}) = 0.023 = B$$

$$SD(\text{total}) = 0.34 = C \quad (7.3 \% \text{ CV, 0-60 days})$$

$$A^2 + B^2 = C^2 \text{ or,}$$

$$SD^2_{\text{sample}} + SD^2_{\text{instrument}} = SD^2_{\text{total}}$$

Today:

$$SD^2_{\text{sample}} + (0.023)^2_{\text{instrument}} = (0.34)^2_{\text{total}},$$


or

$$SD^2_{\text{sample}} = 0.339$$

On-farm, if CV% is 1.5%:

$$(x)^2_{\text{sample}} + (0.07)^2_{\text{instrument}} = (0.34)^2_{\text{total}}, \text{ or}$$

$$SD^2_{\text{sample}} = 0.33$$



lattec

Conclusion:

- Total variation around 7-10 % at the beginning of lactation seems to be the case.
- A number of error sources are abolished or minimized with on-farm analyses, a few increases.
- We should use all information available to improve the total setup, - in stead of focusing only on parts of the system
- Work to establish technical specifications for on-farm analyzers is needed


lattec

Acknowledgements:

- Dr Dave Kelton, Dr Karen Hand, University of Guelph, Canada
- Dr Carsten Ridder, Lattec, Denmark
- Dr Peter Løvendahl et al., Danish Institute of Animal Science

lattec

Standard deviation (prot) over 5 days



lattec

lattec

Focus is on start of lactation!

| cow | SD(5) | | CV%(5) | av fat% |
|----------|---------|---------|--------|---------|
| | 60 days | 60 days | | |
| 244-4 | 0.37 | 6.10 | 5.68 | |
| 485-1 | 0.27 | 5.33 | 3.79 | |
| 392-1 | 0.27 | 6.95 | 3.94 | |
| 281-5 | 0.33 | 6.50 | 4.75 | |
| 480-1 | 0.25 | 6.05 | 3.43 | |
| 469-1 | 0.41 | 9.88 | 4.02 | |
| 8202-5 | 0.41 | 7.50 | 5.23 | |
| 499-1 | 0.28 | 6.13 | 4.54 | |
| 240-3 | 0.50 | 8.16 | 5.36 | |
| 287-5 | 0.29 | 5.11 | 4.44 | |
| 809-1 | 0.27 | 5.58 | 5.35 | |
| 390-Prot | 0.14 | 4.92 | 3.49 | |
| | 0.34 | | 4.50 | |

For all lactation:

Average CV (DK):
5.25, (or SD = 0.24)

Average CV (CND):
6.40, (or SD = 0.29)

Average CV Parlor1:
6.43

Average CV Parlor 2:
6.54

lattec